### Description

The *TorEx Lung Perfusion Kit* contains single-use components required during ex vivo lung perfusion using the *TorEx Lung Perfusion System* (906-0004). The kit is provided sterile and includes:

- *TorEx Lung Perfusion Organ Chamber* (902-0003)
- *TorEx Lung Perfusion LA Cannulae*, provided in three sizes (903-0006, 903-0007, 903-0008)
- *TorEx Lung Perfusion PA Cannula* (903-0009)
- *TorEx Lung Perfusion Rescue Cannula Kit* (903-0010)
  - *TorEx Lung Perfusion LA Rescue Cannula* (903-0011)
  - *TorEx Lung Perfusion PA Rescue Cannula* (903-0012)
- *TorEx Lung Perfusion Drain Kit* (907-0001)
- Ventilation hose (201-0003)
- Ventilation calibration adapter (402-0056)
- Drip tray (500-0001)
- Biohazard trash liner (805-0034)
- Hydrostatic pressure sensing lines (301-0065)
- Uline Economy Cutter (703-0029)
- Cleaning brush (703-0054)
- Zip ties (700-0010)

The *TorEx Lung Perfusion Organ Chamber* holds the donor's lungs in a sterile chamber. It also houses the entire perfusion fluid circuit, including the pump head and the gas exchanger, as well as ports to add, drain and sample the perfusate fluid. It also includes ventilation tubes, filters and sensors.

The *TorEx Lung Perfusion Cannulae* are used to connect the perfusion fluid circuit to the pulmonary artery and to the left atrium cuff of the lungs. There is a PA cannula (for the pulmonary artery) and a LA cannula (left atrium). The LA cannula is provided in multiple sizes. Rescue cannulae are also provided, to be used when the LA or PA cannulae cannot be used because insufficient tissue is left on the LA cuff or pulmonary artery. Every cannula includes a connection for a hydrostatic pressure sensor.

The *TorEx Lung Perfusion Drain Kit* is used to drain used perfusate from the waste reservoir of the *TorEx Lung Perfusion Organ Chamber*.

### Intended Use





The *TorEx Lung Perfusion EVLP Kit* is intended to be used as part of the *TorEx Lung Perfusion System* for temporary continuous closed system warm perfusion of excised donor lungs, during which time the function of the lungs can be assessed for transplantation.

### Warnings and Cautions

- Do not use if the packaging of the *TorEx Lung Perfusion Kit* is damaged.
- Components of the *TorEx Lung Perfusion Kit* should be used as soon as the package is open.
- Components of the *TorEx Lung Perfusion Kit* should be disposed of as biohazardous waste once they are used.

- Components of the *TorEx Lung Perfusion Kit* cannot be resterilized after they have been removed from their sterile package.

### Symbols on labels

Symbol	Meaning
	Type BF applied part
	Refer to Instructions for Use
	Manufacturer
	Date and country of manufacture
	Catalog number (part number)
	Lot number
	Serial number
	Use by date
	Medical device contains biological tissue, cells, or their derivatives, of animal origin
	Do not use if package is damaged
	Temperature limits for storage
	Humidity limits for storage
	Do not reuse
	Single sterile barrier system
	Do not resterilize
	Indicates the item is a medical device
	Caution : Federal (US) law restricts this device to sale by or on the order of a physician.

### Instructions

Please follow the Instructions for Use for the *TorEx Lung Perfusion System* (906-0004-UM-01).

Traferox Technologies Inc.  
3505 Laird Rd. Unit 16  
Mississauga, ON, Canada, L5L 5Y7

Customer Support: +1 905-503-5419  
Email: support@traferox.com  
www.traferox.com

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