en: Instructions for Use

TorEx Lung Perfusate

For assessment of isolated lungs

Product Description

The TorEx Lung Perfusate is a physiological salt solution containing (70 g/L) human serum albumin and (5 g/L) Dextran 40. The electrolyte composition, pH and human serum albumin simulate key properties of human blood plasma. Dextran 40 counteracts tissue oedema and protects the microvasculature against post-ischaemic reperfusion injury. Human serum albumin provides oncotic pressure preventing oedema.

The solution is sterile (aseptic processing) and intended for single use only.

The bottle is made of PET-G and equipped with a PE screw cap lined with a silicone septum closure which facilitates aseptic transfer of the solution. The outer carton box is closed with a tamper-evident sticker.

TorEx Lung Perfusate enables safe Ex Vivo Lung Perfusion (EVLV).

The intended patient population is adult patients in need of a lung transplantation.

Clinical Benefit

TorEx Lung Perfusate enables safe Ex Vivo Lung Perfusion (EVLV) and assessment of isolated lungs, increasing the number of donor lungs available for transplant.

Intended Purpose

The TorEx Lung Perfusate is intended for assessment of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient.

Contraindication

There are no known contraindications.

Precautions

- Do not use TorEx Lung Perfusate if the solution is unclear, the bottle is damaged, the tamper-evident sticker has been broken or if the expiry date has passed.
- TorEx Lung Perfusate solution is intended for single use only and **MUST NOT BE REUSED**. Reuse of TorEx Lung Perfusate solution is not allowed due to the risk of infection.
- TorEx Lung Perfusate solution is not intended for any other use than that stated in the section Intended Purpose.
- TorEx Lung Perfusate solution is not intended for infusion.
- All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for HBsAg, HIV-1/HIV-2, HAV, HBV, HCV and Parvo B19.
- No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Treatment of children: There is no data on the use of TorEx Lung Perfusate on children.

Pregnancy: There is no data on the use of TorEx Lung Perfusate during pregnancy.

Lactation: It is not known to what extent TorEx Lung Perfusate passes into breast milk

Adverse Effects

When administered systemically, human serum albumin has been associated with rare allergic reactions (<1 in 1000) such as urticaria, fever, chills, pruritus, and anaphylaxis.

Dextran has also been associated with rare anaphylactic reactions (<1 in 1000). However, no such reactions have been reported with either of these substances when used for ex vivo organ or tissue perfusion or preservation.

Any invasive procedure carries a risk of infection.

Storage Conditions

TorEx Lung Perfusate solution shall be stored in the cardboard box at +2°C to +8°C with the bottle in an upright position.

Instructions

Lung evaluation set-up

The responsibility for correct clinical use and technique rests with the user. The Instructions are only provided as a suggestion for procedure, the user must on the basis of her or his medical training and experience evaluate the suitability of this procedure.

Perfusion pressure is closely monitored and should never exceed 20 mmHg. Since the lungs are generally cold before perfusion, initial flow rates must be very low, gradually increasing as the temperature increases. The flow should never exceed 4 L/min.

Functional assessment can begin when the circuit has reached the steady state at 37°C, monitoring arterial and venous blood gases, pulmonary vascular resistance, ventilation compliance and other relevant parameters.

1. Read and understand the instructions below before proceeding.
2. Assemble a sterile extra corporal circulation system comprising pump, oxygenator, heater/cooler, fluid reservoir and leucocyte filter. A ventilator is required to ventilate the lung(s). Keep a blood gas analyzer in close proximity to the system, to allow for repeated measurement of blood gases, pH and other relevant parameters.
3. Remove the bottle from the carton and check that:
 - the tamper-evident sticker on the carton box is intact,
 - the bottle is not damaged or leaking, and
 - the expiry date has not passed.

Perfusate should be dispensed as required by the ex vivo perfusion technique.

4. Disinfect the septum closure of the screw cap and transfer the solution aseptically to the fluid reservoir in the extra corporal circulation system. Circulate the solution and de-air the system.
5. If the system is not primed properly, repeat with a new bottle of TorEx Lung Perfusate from procedure 3.
6. Acellular TorEx Lung Perfusate can be used as-is for lung perfusion.
7. Circulate the final solution until appropriate baseline values of pH, pCO₂, pO₂, and any other relevant parameters have been established. Ensure that the fluid reservoir always contains sufficient amount of fluid in order to prevent air entering the EVLV circuit.
8. Once appropriate values have been established, proceed with connecting the lungs, flushed and stored in low potassium Dextran preservation solution, to the extra corporal circulation system and progressively increase perfusate temperature to normothermia before commencing evaluation of lung function according to established clinical routines. Do not commence ventilation of the lung(s) until a temperature of 32°C has been reached in the perfusate exiting the lung(s). Clinical judgement must in each case be the decisive factor for when to end EVLV and if to transplant or not.
9. Dispose of the bottles and any residual solution after the procedure according to the hospital's standard routines. Note that the solution contains human serum albumin.

Perfusate should be dispensed as required by the ex vivo perfusion technique.

6. Acellular TorEx Lung Perfusate can be used as-is for lung perfusion.
7. Circulate the final solution until appropriate baseline values of pH, pCO₂, pO₂, and any other relevant parameters have been established. Ensure that the fluid reservoir always contains sufficient amount of fluid in order to prevent air entering the EVLV circuit.

Once appropriate values have been established, proceed with connecting the lungs, flushed and stored in low potassium Dextran preservation solution, to the extra corporal circulation system and progressively increase perfusate temperature to normothermia before commencing evaluation of lung function according to established clinical routines. Do not commence ventilation of the lung(s) until a temperature of 32°C has been reached in the perfusate exiting the lung(s). Clinical judgement must in each case be the decisive factor for when to end EVLV and if to transplant or not.

Dispose of the bottles and any residual solution after the procedure according to the hospital's standard routines. Note that the solution contains human serum albumin.

Additives

Additives may be added to the reservoir at the determination of the clinician. Those additives can be broad spectrum antibiotics, corticosteroids such as methylprednisolone, heparin and when indicated antifungal medication. The dosage depends on the substance(s) used and the decision of the clinician. Be careful to supplement with the additives also in the replacement solution if part of TorEx Lung Perfusate is exchanged during perfusion. Additional additives have not been thoroughly investigated but other additives, mixable with plasma, might be used based on the experience of the clinician.

Symbols

	Refer to Instructions for Use		Do not use if package is damaged
	Manufacturer		Keep away from sunlight
	Catalog number (part number)		Contains human blood or plasma derivatives
	Lot number		Contains medicinal substance
	Date and country of manufacture		Contains Dextran 40
	Use by date		Temperature limitation
	Do not reuse		Medical Device
	Do not re-sterilize		Store in upright position
	Single sterile barrier system with protective packaging outside	Rx only	Caution: Federal US law restricts this device to sale by or on the order of a physician.
	Sterilized using aseptic processing techniques		

Contact Information:

Traferox Technologies Inc.
3505 Laird Rd. Unit 16
Mississauga, ON, Canada, L5L 5Y7
Canada

Technical Support
Phone: +1 905-503-5419 option 1
Email: support@traferox.com

TorEx® is a registered trademark of Traferox Technologies Inc.

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TorEx® est une marque de commerce enregistrée de Traferox Technologies Inc.